

Vibrio tubiashii

区分	手法名 (文献)	プライマー		反応温度条件	増幅産物 bp	備考	推奨度
		名称	配列(5'-3')				
PCR	VS3/VS4 (Sathyamoorthy et al. 2011)	VS3	GCTCGAAGCTCTTTGGATTG	94°C5分→(94°C30秒、55°C1分、72°C1分)×26サイクル→72°C7分	250	—	—
		VS4	TCCGCCATATCCCTACTCAC				

文献

Sathyamoorthy, V., Datta, A.R., Lee, C.J., Kothary, M.H., McCardell, B.A. and Tall, B.D. (2011). Cloning and partial characterization of a novel hemolysin gene of *Vibrio tubiashii* and the development of a PCR-based detection assay. *Can. J. Microbiol.* 57, 714-721.